



3 - Clamp Head Mask
3HM07-00



4 - Clamp Neck & Shoulder Mask
4HM08-00



5 - Clamp Head, Neck & Shoulder Mask
5HM09-00



Prone Breast Mask
PBM-33



Breast Mask
BM10-00



4 - Clamp Pelvic Mask
4PM11-00



6 - Clamp Pelvic Mask
6PM12-00

Thermoplastic Masks



Immobilization devices are used in radiation therapy to restrict patient movement during treatment. By restricting movement, the immobilization device helps to ensure that the radiation beam is precisely targeted to the intended treatment area while minimizing radiation exposure to healthy tissues. Panacea's skin-safe IMPACT masks are non-sticky thermoplastics with uniform perforation. The masks can be moulded on patient's body with ease with instant hardening and lesser shrinkage. IMPACT masks exert less pressure on the patient's skin and thus enable better patient comfort. It also safeguards the reproducibility and stability to deliver high precision patient immobilization improving the overall treatment experience.

Salient features:

- ✓ Uniform perforation
- ✓ Maximum stability
- ✓ Proficient daily reproducibility
- ✓ No hostile effects on skin
- ✓ Available optimum thicknesses:
2.4 mm and 3.2 mm